

	STANDARD SUPASCOOTA	HEAVY DUTY SUPASCOOTA
Wheel size:	200 mm	200 mm
Ground clearance:	73 mm	70 mm
Length:	930 mm	945 mm
Width:	535 mm	565 mm
Weight (without seat):	13 kg	18.3 kg
Weight of seat:	3.2 kg	7.2 kg
Max loading:	80 kg	130 kg
Turn radius	780 mm	780 mm
Motor output	DC 12v - 180 Watt	2 - 180 Watt
Transmission:	Worm drive gearbox	Worm drive gearbox
Reverse:	Forward/reverse switch	Forward/reverse switch
Brakes:	Dynamic braking and electromagnetic brake	Dynamic braking and electromagnetic brake
Battery (SLA):	12V - 20 ah (6.6 kg)	24V - 15 ah (9.2 kg)
Charger:	Off-board 12 volt 3 amp; (Optional in-car charging lead available)	24 volt 2 amp
Lithium battery (optional):	14.4V / 20ah 2.3kgs (4-5 times life cycle of SLA battery) Lithium charger included (optional)	-
Max speed SLA battery:	6.5 kph	6.5 kph
Max speed Lithium battery:	7.7 kph	-
Slope grade ability:	Max. recommended 9 degrees (80kg test)	9 degrees with 120kg load 12 degrees with 100kg load
Travel distance:	Up to 10km with standard SLA battery Up to 20km on optional Lithium battery	Up to 16km with standard SLA battery
Arm rests (optional):	Padded, flip back, removable; Distance between armrests 45cm or 50cm	Padded, flip back, removable; Distance between armrests 45cm or 50cm